

# MNova Installation & Activation Guide

Indian Institute of Science, Bangalore



## STEP 1: Test Connection

MNova License service is installed on **server 10.16.4.105**. Ping this address from your local machine (computer or laptop), where you want to install MNova. If it's a successful ping, you can proceed with step 2; otherwise, contact the SERC team. Contact details are provided at the end of this document.

🔊 ping 10.16.4.105

"Alternatively, you may proceed without using the ping command. Simply enter the IP address 10.16.4.105 whenever you are prompted to provide the server name or IP address."

## STEP 2: Download the latest version of MNova.

<https://mestrelab.com/download/mnova/>

## Download the latest Mnova Version



### Latest Manual:

Please note that manuals are version-specific. If you need a manual for a previous version, please go to the bottom of the page and select the version that you are interested in.

[Download the latest manual \(PDF\)](#)

### Minimum system requirements

Windows 7 or higher

Apple OS X 10.13 or higher.

Pentium 300 MHz, 128MB RAM, Video Adapter Super VGA (800 x 600) with X11, OpenSSL and OpenGL libraries.

### Additional Software:

If you do not have system administrator permissions please:

– [Download our per user installer\(64-bits\)](#)

Choose your distro and make sure you also download and install our sign key.



## STEP 3: Installation

**NOTE:** If you are going to install the program on a Windows system, be sure you have administrator privileges, because the installation procedure will install several files in your system folders. If you don't have such privileges, you won't be able to install or uninstall the program correctly.


- Double click setup file 7 next



- Accept license agreement 7 next



- Install



## STEP 4: Activation

Run the application. The following dialogue, notifying you that some licenses were not found, might appear:



Drag & drop the license file to MNova in order to activate it. Alternatively, you can also follow the Menu →Help→ Get-Install Licenses and click on the 'Install' button.

This will open the ‘Select License File’ dialog (see the picture below). Navigate on your Windows browser until you find the license file (or a .zip file containing several licenses), double-click on it, and select it.



OR

Copy the license file and save it to

**C:\Program Files\Mestrelab Research S.L\MestReNova\licenses**

**For MNova License Key Request:** For your personal license key for the MNova software, kindly email [services.library@iisc.ac.in](mailto:services.library@iisc.ac.in) with the “MNova License Key Request” subject line. The library team will promptly furnish you with the activation details.

**License Server Queries:** The license server for MNova is operational at the SERC facility. Should you encounter any issues or have inquiries regarding the license server, please contact:

Email: [sysadmin.serc@iisc.ac.in](mailto:sysadmin.serc@iisc.ac.in)

Ext: 119